

Rabbit Anti-Mouse Vcam-1/Cd106 Monoclonal Antibody For Fc - P29533

Item Number: CM0005484

Introduction

-conjugated rabbit monoclonal antibody specific for mouse VCAM-1 (CD106), validated for flow cytometry. Targets leukocyte-endothelial adhesion and immune surveillance pathways. Recombinant immunogen covers extracellular domain.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse VCAM-1/CD106 mAb
Remarks/Alias	CD106; Vcam-1
Species	Mouse
Gene ID (Human)	7412
Gene ID	22329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-310 of mouse VCAM-1/CD106 (NP_035823.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P29533
Protein Weight	38kDa/81kDa
Shipping	Ice bag