

Rabbit Anti-Mouse Vista Polymab® For Fc - Q9D659

Item Number: CM0005267

Introduction

-conjugated rabbit monoclonal antibody targeting mouse VISTA (Dies1/PD-1H). Validated for flow cytometry (FC). Useful for studying T-cell suppression, cancer immunotherapy, and autoimmune disease.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse VISTA PolymAb®
Remarks/Alias	Dies1; PD-1H; VISTA; 4632428N05Rik
Species	Mouse
GeneID (Human)	64115
GeneID	74048
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-191 of mouse VISTA (NP_083008.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9D659
Protein Weight	34kDa
Shipping	Ice bag