

Rabbit Anti-Mouse/Human Cd55 Mab For Fc - Q61475

Item Number: CM0007832

Introduction

-conjugated rabbit monoclonal antibody targeting mouse/human CD55 (DAF) for flow cytometry. Recognizes complement decay-accelerating factor, GPI-anchored, 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Human CD55 mAb
Remarks/Alias	Daf; Daf1; Daf-GPI; GPI-DAF
Species	Mouse
Gene ID (Human)	1604
Gene ID	13136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-361 of mouse CD55/DAF (NP_034146.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q61475
Protein Weight	42KDa
Shipping	Ice bag