

Rabbit Anti-Mouse/Human Cd55 Monoclonal Antibody For Flow Cytometry - Q61475

Item Number: CM0009231

Introduction

-conjugated rabbit monoclonal antibody targeting mouse and human CD55/DAF. Validated for flow cytometry, this 42 kDa GPI-anchored complement regulator reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Human CD55 mAb
Remarks/Alias	Daf; Daf1; Daf-GPI; GPI-DAF
Species	Mouse
Gene ID (Human)	1604
Gene ID	13136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-361 of mouse CD55/DAF (NP_034146.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q61475
Protein Weight	42KDa
Shipping	Ice bag