

## Rabbit Anti-Mouse/Rat Cd27 Monoclonal Antibody For Fc - P41272

Item Number: CM0012089

### Introduction

Rabbit anti-Mouse/Rat CD27 monoclonal antibody for flow cytometry. Targets CD27 (TNFRSF7), a costimulatory immune-checkpoint receptor. Validated in mouse and rat. Suitable for T cell, NK cell, and B cell studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse/Rat CD27 mAb                                                                                      |
| Remarks/Alias    | S152; Tp55; Tnfrsf7                                                                                                 |
| Species          | Mouse                                                                                                               |
| GeneID (Human)   | 939                                                                                                                 |
| GeneID           | 21940                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-182 of mouse CD27 (NP_001028298.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | FC                                                                                                                  |
| Cross Reactivity | Mouse, Rat                                                                                                          |
| SWISS            | P41272                                                                                                              |
| Protein Weight   | 28kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |