

Rabbit Anti-Mouse/Rat Cd27 Monoclonal Antibody For Fc - P41272

Item Number: CM0010987

Introduction

-conjugated rabbit monoclonal antibody against mouse/rat CD27, validated for flow cytometry (FC). Recognizes CD27 antigen, a costimulatory immune checkpoint receptor critical for T-cell activation. Molecular weight: 28 kDa. Species reactivity: mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD27 mAb
Remarks/Alias	S152; Tp55; Tnfrsf7
Species	Mouse
GeneID (Human)	939
GeneID	21940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-182 of mouse CD27 (NP_001028298.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P41272
Protein Weight	28kDa
Shipping	Ice bag