

Rabbit Anti-Mouse/Rat Cd31/Pecam1 Monoclonal Antibody For Fc - Q08481

Item Number: CM0008444

Introduction

-conjugated rabbit monoclonal antibody targeting mouse/rat CD31/PECAM1. Validated for flow cytometry (FC), cross-reacts with mouse and rat. Immunogen: amino acids 18-220 of mouse CD31.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse/Rat
GenelD (human)	5175
GenelD	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	81 kDa
Shipping	Ice bag