

Rabbit Anti-Mouse/Rat Cd31/Pecam1 Monoclonal Antibody For Flow Cytometry - Q08481

Item Number: CM0010748

Introduction

-conjugated rabbit monoclonal anti-CD31/PECAM1 antibody for flow cytometry. Detects mouse and rat CD31, a cell adhesion molecule crucial for leukocyte transendothelial migration and endothelial cell-cell interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
Gene ID (Human)	5175
Gene ID (Mouse)	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	69kDa/77kDa/80kDa/81kDa
Shipping	Ice bag