

Rabbit Anti-Mouse/Rat Cd54 Mab - P13597

Item Number: CM0012803

Introduction

-conjugated rabbit monoclonal antibody targeting mouse/rat CD54 (ICAM-1) for flow cytometry. Detects a 59 kDa cell adhesion molecule critical for leukocyte migration and immune response.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD54 mAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2
Species	Mouse
Gene ID (Human)	3383
Gene ID	15894
Immunogen	Recombinant Protein / Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse CD54(NP_034623.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P13597
Protein Weight	59kDa
Shipping	Ice bag