

Rabbit Anti-Mouse/Rat Cd54 Monoclonal Antibody For Flow Cytometry - P13597

Item Number: CM0008395

Introduction

-conjugated rabbit anti-mouse/rat CD54 (ICAM-1) monoclonal antibody for flow cytometry. Detects ~59 kDa protein, cross-reacts with Mouse and Rat. Ideal for studying leukocyte adhesion and trans-endothelial migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD54 mAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2
Species	Mouse
Gene ID (Human)	3383
Gene ID (Mouse)	15894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse CD54 (NP_034623.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P13597
Protein Weight	59kDa
Shipping	Ice bag