

Rabbit Anti-Myh7/β-Mhc Polyclonal Antibody For Wb, If/Icc, Elisa - P12883

Item Number: CM0030330



Introduction

Rabbit polyclonal antibody targeting MYH7/β-MHC (Myosin-7). Validated for WB, IF/ICC, and ELISA applications. Reacts with human, mouse, and rat samples. Ideal for cardiac and skeletal muscle research, including cardiomyopathy studies.

[Learn More](#)

Parameter	Value
Product Name	MYH7/β-MHC Rabbit pAb
Remarks/Alias	CMH1; MPD1; SPMD; SPMM; CMD1S; MYHCB; CMYP7A; CMYP7B; MYH7/β-MHC
Species	Human
GeneID (Human)	4625
GeneID	4625
Immunogen	Synthetic peptide / A synthetic peptide corresponding to a sequence within amino acids 1700-1800 of human MYH7/β-MHC (NP_000248.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12883
Protein Weight	223kDa
Shipping	Ice bag