

Rabbit Anti-Nhp Igd Mab - Anti-Nhp Igd

Item Number: CM0010799

Introduction

Rabbit anti-NHP IgD mAb is a monoclonal antibody for flow cytometry, reactive with cynomolgus and rhesus monkey IgD. Ideal for non-human primate immunophenotyping research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag