

Rabbit Anti-Nhp Igd Mab - Nhp Igd Monoclonal Antibody

Item Number: CM0003985

Introduction

Rabbit anti-NHP IgD monoclonal antibody for flow cytometry. Designed for cynomolgus and rhesus monkey immunoglobulin D detection in non-human primate research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag