

Rabbit Anti-Nhp IgD Monoclonal Antibody For Flow Cytometry - Nhp IgD

Item Number: CM0003972

Introduction

-conjugated rabbit monoclonal antibody targeting non-human primate IgD, suitable for flow cytometry (FC) with cross-reactivity to cynomolgus and rhesus monkey samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag