


Rabbit Anti-Ogfod1 Polyclonal Antibody For Wb, Elisa - Q8N543

Item Number: CM0018956



Introduction

Rabbit polyclonal antibody against human OGFOD1 (TPA1), validated for Western blot and ELISA. Targets the 63kDa prolyl 3-hydroxylase involved in translation termination and stress granule formation.

[Learn More](#)

Parameter	Value
Product Name	OGFOD1 Rabbit pAb
Remarks/Alias	TPA1; OGFOD1
Species	Human
GeneID (Human)	55239
GeneID	55239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human OGFOD1 (NP_060703.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N543
Protein Weight	63kDa
Shipping	Ice bag