


Rabbit Anti-Orc3 Polyclonal Antibody For Wb, Elisa - Q9Ubd5

Item Number: CM0016776



Introduction

High-quality rabbit polyclonal antibody targeting human ORC3, validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Suitable for investigating origin recognition complex function, DNA replication initiation, and epigenetic regulation.

[Learn More](#)

Parameter	Value
Product Name	ORC3 Rabbit pAb
Remarks/Alias	LAT; ORC3L; LATHEO; ORC3
Species	Human
Gene ID (Human)	23595
Gene ID	23595
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ORC3 (NP_036513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBD5
Protein Weight	82kDa
Shipping	Ice bag