

Rabbit Anti-P70 S6 Kinase 2 Monoclonal Antibody For Wb - Q9Ubs0

Item Number: CM0006227

Introduction

This rabbit monoclonal antibody detects human, mouse, and rat ribosomal protein S6 kinase beta-2 (p70 S6K2), a key mTOR pathway kinase. Validated for WB, IHC-P, IF/ICC, and ELISA.

[Learn More](#)

Parameter	Value
Product Name	p70 S6 Kinase 2 Rabbit mAb
Remarks/Alias	KLS, SRK, S6K2, S6KB, STK14B, S6KI(2), S6Kbeta, p70S6Kb, P70-beta, S6K-beta2, P70-beta-1, P70-beta-2, p70(S6K)-beta, p70 S6 Kinase 2
Species	Human
Gene ID (Human)	6199
Gene ID	6199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human p70 S6 Kinase 2 (Q9UBS0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBS0
Protein Weight	53 kDa
Shipping	Ice bag