

Rabbit Anti-Perk Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nzj5

Item Number: CM0019889



Introduction

Rabbit polyclonal antibody targeting human PERK (EIF2AK3/PEK), suited for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PERK. A key tool for unfolded protein response and stress signaling research.

[Learn More](#)

Parameter	Value
Product Name	PERK Rabbit pAb
Remarks/Alias	PEK; WRS; PERK
Species	Human
GeneID (Human)	9451
GeneID	9451
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 867-1116 of human PERK (NP_004827.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ5
Protein Weight	125kDa
Shipping	Ice bag