

Rabbit Anti-Pig Cd138 Monoclonal Antibody For Fc - K7Gnf5

Item Number: CM0009367

Introduction

-conjugated rabbit monoclonal antibody against pig CD138 (syndecan-1, 32 kDa proteoglycan) for flow cytometry. Targets extracellular domain (aa 23-250, UniProt K7GNF5).

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Pig CD138 mAb
Remarks/Alias	SDC1
Species	Pig
Gene ID (Human)	6382
Gene ID	100519770
Immunogen	Recombinant fusion protein containing amino acids 23-250 of pig CD138 (NP_001230119.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
SWISS	K7GNF5
Protein Weight	32kDa
Shipping	Ice bag