

## Rabbit Anti-Pig Cd27 Monoclonal Antibody For Flow Cytometry - L7Qhb1

Item Number: CM0005795

### Introduction

-conjugated rabbit monoclonal antibody against pig CD27 for flow cytometry. Recognizes pig CD27, a costimulatory immune-checkpoint receptor. Ideal for T-cell activation, memory, and B-cell surveillance studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Pig CD27 mAb                                                                                          |
| Species          | Pig                                                                                                               |
| GeneID (Human)   | 939                                                                                                               |
| GeneID           | 100520023                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-180 of pig CD27 (NP_001265692.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Pig                                                                                                               |
| SWISS            | L7QHB1                                                                                                            |
| Protein Weight   | 27kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |