

Rabbit Anti-Pig Cd3G Monoclonal Antibody - Cd3G

Item Number: CM0011785

Introduction

-conjugated rabbit monoclonal antibody targeting pig CD3G for flow cytometry applications. Enables reliable detection of CD3G in porcine samples, suitable for T-cell immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD3G Monoclonal Antibody
Species	Pig
Gene ID (Human)	917
Gene ID	494013
Immunogen	Recombinant protein of pig CD3G
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
Shipping	Ice bag