

Rabbit Anti-Pig Cd3G Monoclonal Antibody - Cd3G

Item Number: CM0011730

Introduction

-conjugated rabbit monoclonal antibody targeting pig CD3G for flow cytometry. Specific to porcine samples, this antibody aids T-cell surface marker analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD3G mAb
Species	Pig
Gene ID (Human)	917
Gene ID	494013
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-120 of pig CD3G. (NP_001008686.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
Shipping	Ice bag