

Rabbit Anti-Pig Cd4 Mab - Cd4

Item Number: CM0008585

Introduction

Rabbit anti-Pig CD4 monoclonal antibody for flow cytometry. Validated for pig samples, immunogen from pig CD4 extracellular domain. Reliable tool for T-cell research in swine models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Pig
Gene ID (Human)	920
Gene ID	404704
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
Shipping	Ice bag