

Rabbit Anti-Pig Cd4 Monoclonal Antibody For Flow Cytometry - Cd4

Item Number: CM0007981

Introduction

-conjugated rabbit monoclonal antibody targeting pig CD4, designed for flow cytometry applications. Specifically reacts with porcine CD4, supporting immune cell profiling in swine research models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Pig
Gene ID (Human)	920
Gene ID	404704
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-395 of pig CD4(NP_001001908.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
Shipping	Ice bag