

Rabbit Anti-Pig Cd8A Monoclonal Antibody For Flow Cytometry - A0A287A9Y5

Item Number: CM0010808

Introduction

-conjugated rabbit anti-pig CD8a monoclonal antibody specifically designed for flow cytometry. Targets the T-cell surface glycoprotein CD8a (SWISS: A0A287A9Y5) with a predicted molecular weight of 26kDa. Suitable for detecting cytotoxic T cells in porcine samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD8a mAb
Species	Pig
Gene ID (Human)	925
Gene ID	396627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-184 of pig CD8a (XP_005662451.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
SWISS	A0A287A9Y5
Protein Weight	26kDa
Shipping	Ice bag