

Rabbit Anti-Ppbp Polyclonal Antibody For Wb And Elisa - P02775

Item Number: CM0019852



Introduction

Rabbit polyclonal antibody against human PPBP, validated for Western blot and ELISA. Reacts with human PPBP (CXCL7), a chemokine involved in neutrophil activation, antimicrobial defense, and growth factor activities. Suitable for cytokine research.

[Learn More](#)

Parameter	Value
Product Name	PPBP Rabbit pAb
Remarks/Alias	PPBP; B-TG1; Beta-TG; CTAP-III; CTAP3; CTAPIII; CXCL7; LA-PF4; LDGF; MDGF; NAP-2; PBP; SCYB7; TC1; TC2; TGB; TGB1; THBGB; THBGB1; platelet basic protein
Species	Human
Gene ID	5473 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-128 of human PPBP (NP_002695.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02775
Protein Weight	13kDa
Shipping	Ice bag