

## Rabbit Anti-Puromycin Mab - Puromycin

Item Number: CM0012423

### Introduction

-conjugated rabbit monoclonal anti-Puromycin antibody for flow cytometry (intracellular). Species-independent detection for puromycin labeling assays. Suitable for IVD raw material use.

[Learn More](#)

| Parameter        | Value                                         |
|------------------|-----------------------------------------------|
| Product Name     | Rabbit anti-puromycin mAb                     |
| Remarks/Alias    | Puromycin                                     |
| Gene ID (Human)  | CAS:58-58-2                                   |
| Gene ID          | CAS:58-58-2                                   |
| Immunogen        | Chemical compounds corresponding to puromycin |
| Source           | Rabbit                                        |
| Category         | Monoclonal Antibodies                         |
| Application      | FC (intra)                                    |
| Cross Reactivity | Species independent                           |
| Shipping         | Ice bag                                       |