

Rabbit Anti-Puromycin Monoclonal Antibody - Puromycin

Item Number: CM0010539

Introduction

-conjugated rabbit monoclonal antibody against puromycin for intracellular flow cytometry. Species-independent reactivity enables detection across all models. Ideal for protein synthesis monitoring and puromycin incorporation assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-puromycin mAb
Remarks/Alias	Puromycin
Gene ID	CAS:58-58-2
Immunogen	Chemical compounds corresponding to puromycin
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag