

Rabbit Anti-Puromycin Monoclonal Antibody - Puromycin

Item Number: CM0005708

Introduction

High-quality Rabbit anti-puromycin monoclonal antibody for flow cytometry (intracellular). Species-independent detection of puromycin. Ideal for protein synthesis research and puromycin incorporation assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-puromycin mAb
Remarks/Alias	Puromycin
Species	Rabbit
Gene ID	CAS:58-58-2
Immunogen	Chemical compounds corresponding to puromycin
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag