

Rabbit Anti-Ralb Polyclonal Antibody For Wb, Ihc-P, Elisa - P11234

Item Number: CM0028205



Introduction

Rabbit polyclonal antibody against RALB, validated for WB, IHC-P, ELISA. Detects human, mouse, rat RALB. For studies of apoptosis, cell migration, and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	RALB Rabbit pAb
Remarks/Alias	RALB
Species	Human
GeneID (Human)	5899
GeneID	5899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALB (NP_002872.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11234
Protein Weight	23kDa
Shipping	Ice bag