

Rabbit Anti-Rat Cd31 Monoclonal Antibody For Flow Cytometry - Q3Swt0

Item Number: CM0011228

Introduction

-conjugated rabbit monoclonal antibody against rat CD31, validated for flow cytometry. Recognizes a 76 kDa platelet endothelial cell adhesion molecule involved in leukocyte transendothelial migration and cell adhesion.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Rat CD31 mAb
Remarks/Alias	CD31; Pecam; CD31/PECAM1
Species	Rat
GenelD (human)	5175
GenelD	29583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-585 of rat CD31 (NP_113779.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	Q3SWT0
Protein Weight	76kDa
Shipping	Ice bag