


Rabbit Anti-Rat Cd4 Monoclonal Antibody For Flow Cytometry - P05540

Item Number: CM0011475



Introduction

-conjugated rabbit monoclonal antibody against rat CD4 (clone Ly49/Klra1) validated for flow cytometry. Ideal for studying T-cell biology, immune response, and MHC class II interactions. Suitable for rat samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat CD4 mAb
Remarks/Alias	Ly49; Klra1; Ly49s2
Species	Rat
Gene ID (Human)	920
Gene ID	24933
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-394 of rat CD4 (NP_036837.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	P05540
Protein Weight	51 kDa
Shipping	Ice bag