

Rabbit Anti-Rat Icam-1/Cd54 Monoclonal Antibody For Fc - Q00238

Item Number: CM0011579

Introduction

-conjugated rabbit monoclonal antibody targeting rat ICAM-1/CD54 (SWISS: Q00238). Validated for flow cytometry (FC) applications, ideal for studying leukocyte adhesion and trans-endothelial migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat ICAM-1/CD54 mAb
Remarks/Alias	CD54; ICAM; Icam1
Species	Rat
GeneID (Human)	3383
GeneID	25464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-493 of rat ICAM-1/CD54 (NP_037099.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	Q00238
Protein Weight	60kDa
Shipping	Ice bag