

Rabbit Anti-Rat Icam-1/Cd54 Monoclonal Antibody For Flow Cytometry - Q00238

Item Number: CM0004760

Introduction

-conjugated rabbit anti-rat ICAM-1/CD54 monoclonal antibody for flow cytometry applications. This antibody specifically detects rat CD54, a cell adhesion molecule critical in leukocyte-endothelial interactions and immune responses. Suitable for rat models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat ICAM-1/CD54 mAb
Remarks/Alias	CD54; ICAM; Icam1
Species	Rat
GeneID (Human)	3383
GeneID	25464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-493 of rat ICAM-1/CD54 (NP_037099.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	Q00238
Protein Weight	60kDa
Shipping	Ice bag