

Rabbit Anti-Rat Ifn-Gamma Monoclonal Antibody () For Flow Cytometry - P01581

Item Number: CM0008659

Introduction

Rabbit monoclonal antibody against rat IFN-gamma, conjugated to for intracellular flow cytometry (FC). Validated for rat cross-reactivity, suitable for immune response and cytokine research studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat IFN-gamma mAb
Alias	If2f, IFNG2, IFN gamma
Species	Rat
Gene ID (Human)	3458
Gene ID	25712
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-156 of rat IFN gamma (NP_620235.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Rat
SWISS	P01581
Protein Weight	18 kDa
Shipping Condition	Ice bag