


Rabbit Anti-Rps6Ka6 Polyclonal Antibody For Wb, Elisa - Q9Uk32

Item Number: CM0015674



Introduction

Rabbit polyclonal antibody targeting RPS6KA6 (RSK4) for WB and ELISA. This constitutively active serine/threonine kinase participates in p53-dependent cell growth arrest, suitable for human and mouse research.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA6 Rabbit pAb
Remarks/Alias	RSK4, RSK-4, p90RSK6, PP90RSK4, S6K-alpha-6, RPS6KA6
Species	Human
Gene ID	27330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 646-745 of human RPS6KA6 (NP_055311.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UK32
Protein Weight	84kDa
Shipping	Ice bag