

Rabbit Anti-Shp/Nr0B2 Polyclonal Antibody For Wb And Elisa - Q15466

Item Number: CM0024298



Introduction

Rabbit polyclonal antibody against human SHP/NR0B2 (SWISS: Q15466). Validated for WB and ELISA, cross-reactive with mouse and rat. Ideal for nuclear receptor signaling and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SHP/NR0B2 Rabbit pAb
Remarks/Alias	SHP; SHP1; SHP/NR0B2
Species	Human
Gene ID (Human)	8431
Gene ID	8431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-124 of human SHP/NR0B2 (NP_068804.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse,Rat
SWISS	Q15466
Protein Weight	28kDa
Shipping	Ice bag