

Rabbit Anti-Thumpd3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bv44

Item Number: CM0015821

Introduction

Rabbit polyclonal antibody targeting THUMPD3, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for tRNA modification and methylation studies.

[Learn More](#)

Parameter	Value
Product Name	THUMPD3 Rabbit pAb
Remarks/Alias	THUMPD3
Species	Human
Gene ID (Human)	25917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human THUMPD3 (NP_056268.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
UniProt Accession	Q9BV44
Protein Weight	57kDa
Shipping	Ice bag