

Rabbit Anti-Zic1 Polyclonal Antibody For Wb, Elisa - Q15915

Item Number: CM0021880

Introduction

Rabbit polyclonal antibody targeting ZIC1, suitable for WB and ELISA. Reacts with human, mouse, rat. ZIC1 is a zinc finger transcription factor crucial for neurogenesis and CNS development.

[Learn More](#)

Parameter	Value
Product Name	ZIC1 Rabbit pAb
Remarks/Alias	ZIC; CRS6; BAIDCS; ZNF201; ZIC1
Species	Human
Gene ID	7545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-210 of human ZIC1 (NP_003403.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15915
Protein Weight	48kDa
Shipping	Ice bag