

Rabbit Anti-Znf217 Polyclonal Antibody For Wb, If/Icc, Elisa - 075362

Item Number: CM0026999



Introduction

Rabbit polyclonal antibody against human ZNF217 for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. ZNF217 is a zinc finger transcription repressor that promotes cell proliferation and AKT phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	ZNF217 Rabbit pAb
Remarks/Alias	ZABC1; ZNF217
Species	Human
Gene ID (Human)	7764
Gene ID	7764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ZNF217 (NP_006517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075362
Protein Weight	115kDa
Shipping	Ice bag