

Rabbit Igg Isotype Control - Monoclonal Rabbit Igg Isotype Control

Item Number: CM0008897

Introduction

This monoclonal rabbit IgG isotype control conjugated to is designed for flow cytometry applications, serving as a reliable negative control to distinguish specific from non-specific antibody binding in cellular analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit IgG isotype control
Immunogen	Synthetic peptide, A synthesized peptide derived from rabbit IgG isotype control.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag