

Rabbit Igg Isotype Control - Rabbit Igg Isotype Control

Item Number: CM0008880

Introduction

Rabbit IgG isotype control is a monoclonal antibody sourced from rabbit, designed for flow cytometry (FC) applications. It serves as a negative control to assess nonspecific binding and background staining in immunoassays.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit IgG isotype control
Immunogen	Synthetic Peptide A synthesized peptide derived from rabbit IgG isotype control
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag