

Rabbit Igg Isotype Control -

Item Number: CM0004379

Introduction

Rabbit IgG isotype control is a monoclonal antibody for flow cytometry (FC) applications. It serves as a negative control to differentiate non-specific background signal from specific antibody staining. Suitable for research and diagnostic use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit IgG isotype control
Immunogen	Synthetic Peptide A synthesized peptide derived from rabbit IgG isotype control
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag