

Rabbit Monoclonal Antibody To Sars-Cov-2 Spike S1 For Wb/If/Elisa - Spike Glycoprotein

Item Number: CM0002204

Introduction

High-quality rabbit monoclonal antibody targeting SARS-CoV-2 Spike S1 glycoprotein. Validated in WB, IF/ICC, and ELISA applications, ideal for virology research and diagnostic assay development.

[Learn More](#)

Parameter	Detail
Product Name	SARS-CoV-2 Spike S1 mAb
Remarks/Alias	spike glycoprotein
Species	SARS-COV-2
Gene ID (Human)	SARS-CoV-2
Gene ID	43740568
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
Shipping	Ice bag