

Rabif Rabbit Pab For Wb, If/Icc, Elisa - P47224

Item Number: CM0017642



Introduction

Rabbit polyclonal antibody against RABIF (MSS4, RASGRF3) for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Roles in vesicular transport via guanine nucleotide exchange.

[Learn More](#)

Parameter	Value
Product Name	RABIF Rabbit pAb
Remarks/Alias	MSS4; RASGRF3; RASGRF3; RABIF
Species	Human
Gene ID (Human)	5877
Gene ID	5877
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human RABIF (NP_002862.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47224
Protein Weight	14kDa
Shipping	Ice bag