

Rad51D Rabbit Mab For Wb And Elisa - 075771

Item Number: CM0002111



Introduction

Rabbit monoclonal antibody against human RAD51D (RAD51 homolog D), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for DNA repair and homologous recombination studies.

[Learn More](#)

Parameter	Value
Product Name	Rad51D Rabbit mAb
Remarks/Alias	TRAD; R51H3; BROVCA4; RAD51L3; Rad51D
Species	Human
GeneID (Human)	5892
GeneID	5892
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Rad51D (075771)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075771
Protein Weight	35kDa
Shipping	Ice bag