


Rad52 Rabbit Mab For Wb, Elisa - P43351

Item Number: CM0005259



Introduction

High-quality Anti-Rad52 Rabbit Monoclonal Antibody validated for Western blot and ELISA. Specific to human Rad52, involved in DNA double-strand break repair and recombination. Ideal for DNA damage and repair research.

[Learn More](#)

Parameter	Value
Product Name	Rad52 Rabbit mAb
Remarks/Alias	RAD52; DNA repair protein RAD52 homolog; Rad52
Species	Human
GeneID (Human)	5893
GeneID	5893
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 319-418 of human Rad52 (P43351).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P43351
Protein Weight	13kDa/15kDa/25kDa/46kDa
Shipping	Ice bag