

Rad54L Rabbit Polyclonal Antibody For Wb And Elisa - Q92698

Item Number: CM0020180

Introduction

RAD54L rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for DNA repair and recombination research. Swiss-Prot: Q92698.

[Learn More](#)

Parameter	Value
Product Name	RAD54L Rabbit pAb
Remarks/Alias	HR54; hHR54; RAD54A; hRAD54; RAD54L
Species	Human
GeneID (Human)	8438
GeneID	8438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 648-747 of human RAD54L (Q92698)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92698
Protein Weight	84kDa
Shipping	Ice bag