

RalA Rabbit Monoclonal Antibody For Wb, Elisa - P11233

Item Number: CM0007712



Introduction

Rabbit monoclonal antibody targeting human RalA, validated for Western blot and ELISA. Recognizes human and monkey RalA protein (24 kDa, UniProt P11233), a GTPase involved in cell proliferation and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	RalA Rabbit mAb
Remarks/Alias	RAL; HINCONS; RalA
Species	Human
GeneID (Human)	5898
GeneID	5898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-206 of human RalA (NP_005393.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Monkey
SWISS	P11233
Protein Weight	24 kDa
Shipping	Ice bag