

Rala Rabbit Polyclonal Antibody For Wb And Elisa - P11233

Item Number: CM0015166



Introduction

Rabbit polyclonal antibody against human RALA (Ras-related protein Ral-A), validated for Western blot and ELISA; cross-reacts with mouse and rat. Suitable for GTPase signaling, membrane trafficking, and cell division research. 24 kDa.

[Learn More](#)

Parameter	Value
Product Name	RALA Rabbit pAb
Remarks/Alias	RAL; HINCONS; RALA
Species	Human
GeneID (Human)	5898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALA (NP_005393.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11233
Protein Weight	24kDa
Shipping	Ice bag